

The New Versatile Analytical VP-SEM...



The Analytical SEM for Studying Large, Heavy & Tall Samples

- Large sample up to 300mm in diameter
- Observable area up to 203mm in diameter
- Observation and EDS analysis on a sample up to 110mm tall
- EDS, EBSD and WDS simultaneously available
- Versatile port layout for various analytical applications
- 5-axis computer eucentric motorized stage with $-20/+90^{\circ}$ tilt
- One-button selectable High Vacuum / Low vacuum Modes
- Ecological design by clean Turbo Molecular Pumping

The New Versatile Analytical VP-SEM...



The Analytical SEM for Studying Large, Heavy & Tall Samples

- Large sample up to 300mm in diameter
- Observable area up to 203mm in diameter
- Observation and EDS analysis on a sample up to 110mm tall
- EDS, EBSD and WDS simultaneously available
- Versatile port layout for various analytical applications
- 5-axis computer eucentric motorized stage with $-20/+90^\circ$ tilt
- One-button selectable High Vacuum / Low vacuum Modes
- Ecological design by clean Turbo Molecular Pumping